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<p>1 Product Description</p> <p>CAMLOG® implants are endosseous implants and are available in different diameters, lengths, and shapes. They are placed surgically in the bone of the maxilla and mandible and support the attachment of functional and esthetic oral restorations in partially and fully edentulous patients. Prosthetic treatment follows, using single crowns, bridges, full dentures, which are attached to the CAMLOG® implants by means of appropriate devices. The CAMLOG® Implant System contains surgical, prosthetic, and dental laboratory components and tools. There are basically no preferred indications for the different implant configurations.</p>
<p>CAMLOG® implants are suitable for indications for oral endosseous maxillary and mandibular implants in functional and/or esthetic oral rehabilitation of partially and fully edentulous patients. Suggestions for specific applications of the different implant types and sizes may be found in the CAMLOG Compendium 1 Surgery.</p> <p>For further information, see also the "Instruction Manual for the CAMLOG® Implant System", Art. No. J8000.0001 and J8000.0040.</p>
<p>Plattform Switching</p> <p>The CAMLOG® Prosthetic Components for Plattform Switching contain PS Abutments, which are apically reduced in diameter in the connection area to the implant. The implant-abutment interface is oriented in the implant axis direction through the reduced PS Abutment diameter. CAMLOG® Abutments for Plattform Switching (PS) are laser-marked with "PS" for easier identification.</p>
<p>Warning: All prosthetical components PS for Plattform Switching are only to be used in combination with SCREW-LINE implants with K article numbers (K-Series)! If healing caps PS are used during the healing phase, to avoid soft-tissue injury, only prosthetic components PS can be used for the further prosthetical maintainance incl. impression taking PS!</p>

<p>1.1 Universal Abutment PS</p> <p>The Universal Abutment PS, like the Universal Abutment, is used in the fabrication of double-crowns (cone/telescopic) and customized crown and bridge constructions and can be appropriately modified. The color-coded Universal Abutment PS is supplied with an abutment screw and is available in implant diameters 3,8, 4,3, 5,0, and 6,0 mm.</p>
<p>Indications:</p> <ul style="list-style-type: none">Double-crown technique Customized crown and bridge technique
<p>Contraindications:</p> <ul style="list-style-type: none">Primary abutment splinting Restorations with an angulation correction of more than 20° to the implant axis Cast-on technique Single restorations with attached free-end pontic Restorations with a length relationship to the implant length greater than 1:1,25
<p>Materials:</p> <ul style="list-style-type: none">Titanium alloy (Ti6Al4V)

<p>1.2 Temporary Abutment PS</p> <p>The Temporary Abutment PS may be used for esthetic immediate restorations and long-term provisionals. It can be custom-trimmed and should be hand-tightened, only. To maintain stability, it must not be exposed to heavy loading.</p> <p>The Temporary Abutment PS is supplied with an abutment screw and is available in implant diameters 3,8, 4,3, 5,0, and 6,0 mm.</p>
<p>Indications:</p> <ul style="list-style-type: none">Temporary restorations for up to 6 months wearing time

<p>Contraindications:</p> <ul style="list-style-type: none">Wearing time of more than 6 months Primary abutment splinting Single tooth restorations with attached free end pontic Restorations with a length relationship to the implant length greater than 1:1,25
<p>Materials:</p> <ul style="list-style-type: none">PEEK (polyether ether ketone)

<p>1.3 Esthomic® Abutments PS</p> <p>The color-coded Esthomic® Abutments PS are available in straight and 15° angled versions. All Esthomic® Abutments PS are offered as Type A and Type B.</p> <p>Type A and Type B differ in their cam arrangement, with a 60° offset in Type B. This allows six rotational positions, and an optimal prosthesis axis is more easily achieved. The abutments are designed for the fabrication of crown and bridge structures (cement-retained, individually screw-retained) and are delivered with an abutment screw.</p> <p>The Esthomic® Abutment PS are available in implant diameters 3,8, 4,3, 5,0, 6,0 mm and in gingival heights of 1,5 –2,5 mm. Special color-coded selection abutments are available to assist in the selection of the abutments in the dental laboratory (also applicable to Plattform Switching models).</p>
--

<p>Warning! Selection abutments must not be used on patients!</p>
--

<p>Indications:</p> <ul style="list-style-type: none">Customized crown and bridge technique (cement-retained, individually screw-retained) Selection abutments for abutment selection on the laboratory model
<p>Contraindications:</p> <ul style="list-style-type: none">Primary abutment splinting Angulation corrections of more than 20° Cast-on technique
<p>Materials:</p> <ul style="list-style-type: none">Titanium alloy (Ti6Al4V)

<p>1.4 Accessories (lab analogs, screws, instruments)</p> <p>Lab analogs are available for the Abutments for Plattform Switching in the appropriate implant diameters.</p> <p>The lab screw only serves the technician to attach the abutment to the working cast (lab analog) and must be hand-tightened, only. The abutment screw is designed for the final insertion of the abutment into the implant. For implant diameters 3,8 and 4,3 mm, use abutment screw Art. No. J4005.1601 (lab screw Art. No. J4006.1601, thread M 1.6) and for implant diameters 5,0 and 6,0 mm, use abutment screw Art. No. J4005.2001 (lab screw Art. No. J4006.2001, thread M 2,0).</p>

<p>Warning! Selection abutments must not be used on patients!</p>
--

<p>Indications:</p> <ul style="list-style-type: none">Customized crown and bridge technique (cement-retained, individually screw-retained) Selection abutments for abutment selection on the laboratory model

<p>Contraindications:</p> <ul style="list-style-type: none">Primary abutment splinting Angulation corrections of more than 20° Cast-on technique

<p>Materials:</p> <ul style="list-style-type: none">Abutment screw: titanium alloy (Ti6Al4V) Lab screw (brown anodized): titanium alloy (Ti6Al4V)

<p>1.4 Accessories (lab analogs, screws, instruments)</p> <p>Lab analogs are available for the Abutments for Plattform Switching in the appropriate implant diameters.</p> <p>The lab screw only serves the technician to attach the abutment to the working cast (lab analog) and must be hand-tightened, only. The abutment screw is designed for the final insertion of the abutment into the implant. For implant diameters 3,8 and 4,3 mm, use abutment screw Art. No. J4005.1601 (lab screw Art. No. J4006.1601, thread M 1.6) and for implant diameters 5,0 and 6,0 mm, use abutment screw Art. No. J4005.2001 (lab screw Art. No. J4006.2001, thread M 2,0).</p>

<p>The screws are tightened with screwdriver 0,05" hex, manual/ratchet (J5315.0510/0501/0502) or screwdriver 0,05" hex with ISO shaft (J5315.0503/0504). The applicable torques are described in section 3.3.</p>

<p>Materials:</p> <ul style="list-style-type: none">Abutment screw: titanium alloy (Ti6Al4V) Lab screw (brown anodized): titanium alloy (Ti6Al4V)

<p>Warning! Do not mix lab screws with abutment screws! Do not use the lab screw (brown anodized) on patients!</p>

<p>2 General Safety and Precautionary Information</p> <ul style="list-style-type: none">An improper procedure during use of the CAMLOG® prosthetic components can lead to failures, implant damage, bone loss, or esthetically unsatisfactory results. The products may be used only by dentists, physicians, and surgeons trained in the implant system. Use of these products requires specialized knowledge and experience in implant density. CAMLOG regularly offers courses or technical consultations on the use of its products. Your local CAMLOG representative will gladly provide advice. Our products are sold only to physicians/dentists and dental laboratories or on their prescription. Not all components are available in all countries. Technical consultation regarding our products can be obtained orally, in writing, by electronic media, and/or through demonstrations. Such information reflect the current state of knowledge and technology as known at the time of market introduction of the product. They do not release the user from the obligation to personally check the product and its suitability for the proposed purposes, indications, and procedures. Product processing and use occur outside our control and therefore are the responsibility of the user. Any liability for damages caused in this way is excluded. Use the product only in combination with CAMLOG® implants and CAMLOG® tools. We accept no warranty or replacement responsibilities in the event of third-party components being used. The CAMLOG® prosthetic system is part of a larger treatment concept and must be used only with the pertinent original parts and tools according to the suggestions and instructions for use provided by the manufacturer. All components of the CAMLOG® Implant System are carefully matched to one another. The use of third-party components can affect system function and results in the exclusion of any warranty or replacement services offered by ALITEC GmbH/CAMLOG. Tools and system components are designed for specific implant and prosthetic product lines and implant diameters. Use with a different implant or prosthetic product line or a different diameter can lead to mechanical failure of system components, tissue injury, or esthetically unsatisfactory results. Because of this, some implant and prosthetic product lines have some dedicated components and tools. The Temporary Abutment PS, Universal Abutments PS, and Esthomic® Abutments PS are suitable for platform switching. Observe the color markings when selecting tools for the required implant diameter.
--

- We guarantee the flawless quality of our products as part of our sales and delivery conditions. Specifications on existing patents, trademarks, or other intellectual rights are not legally binding.
- During intraoral use, products must be secured in general against aspiration and swallowing.
- Where indications are listed for a particular product, it should be noted that any indications that are not listed are in fact contraindicated.

<p>3 Use</p> <p>3.1 Prosthetic Treatment</p> <p>Final prosthetic treatment with an implant can only be performed after the soft tissue has healed and is free of inflammation. Radiographic follow-up should be performed after 6–12 weeks of healing and before the start of prosthetic treatment. Silicone or polyether impression materials are suitable for taking subsequent impressions of the implant situation.</p>

<p>During the selection of PS Abutments, careful attention should be given to the diameter and angulation of the implant as well as to the gingival height. Angulation corrections of more than 20° should not be performed. No free-end pontic should be added to a free-standing single implant.</p> <p>WARNING! The contact surfaces between the abutment and the implant must not be sandblasted or trimmed!</p> <p>The following points should be taken care of during fabrication of the prosthetic restoration:</p> <ul style="list-style-type: none">Satisfactory load distribution. Tension-free seating of the restoration on the abutments. Correct occlusion. Only specifically recommended materials should be used in prosthetic restorations. Remove the gingiva former and clean the inside of the implant prior to insertion of the abutment in the patient's mouth. The abutment must be precisely seated in the implant. Make sure no soft tissue is trapped and there are no gaps.

<p>3.2 Abutment Selection</p> <p>The abutment is selected by the clinician according to the individual patient case. Abutment selection is dependent on implant diameter and mucosal thickness. Observe color-coding when selecting abutments. In the case of restorations with cementable superstructures, the abutment shoulder should not lie deeper than 1,5 mm below the mucosa so that excess cement can be easily removed.</p>
--

<p>3.3 Insertion of the Abutment into the Implant</p> <p>Remove the healing cap or temporary abutment and clean the inside of the implant prior to insertion of the PS Abutment into the implant.</p>
--

<p>Torques</p> <p>The PS Abutments are attached definitively to the implant using the following torque values:</p>	
<ul style="list-style-type: none">Temporary Abutment PS with abutment screw	hand-tightened
<ul style="list-style-type: none">Universal Abutment PS with abutment screw Esthomic® Abutments PS with abutment screw	20 Ncm
The abutment screws in the Universal und Esthomic® Abutments PS must be retightened after about 5 minutes with the same torque value!	

<p>A control x-ray is advisable to make sure the abutment is correctly seated on the implant.</p>

<p>4 Cleaning, Disinfection and Sterilization</p> <p>The prosthetic components of the CAMLOG® Implant System are supplied non-sterile unless they are explicitly marked as sterile. They must be used one time only and solely on one patient. They must be cleaned and disinfected before each use on a patient and after each use on a patient, e.g., for shipment to the dental laboratory. We recommend an additional sterilization. Exceptions: plastic parts like impression caps for impression posts and bite registration caps must not be sterilized (only cleaned and disinfected).</p> <p>Permissible steam sterilization procedures are fractionated vacuum procedure (with sufficient product drying) and gravitation procedure. Other sterilization procedures are not allowed.</p> <p>Care should be taken that:</p> <ul style="list-style-type: none">The steam sterilization procedure is validated. The maximum steam sterilization temperature is 138 °C/280 °F. The sterilization time (exposure time at sterilization temperature) is at least 20 minutes at 121 °C/250 °F (gravity procedure as well as fractionated vacuum procedure) or 5 minutes at 134 °C/273 °F (gravity procedure as well as fractionated vacuum procedure).
--

<p>Für further information, see also the "Preparation Instructions for the CAMLOG® Implant System", Art. No. J8000.0032.</p>
<p>Attention: All non-sterile packaged CAMLOG® products must not be sterilized in the CAMLOG original packaging!</p>
<p>5 Explanation of Symbols</p> <p>The symbols comply with DIN EN 980. Please refer to the reverse.</p>
<p>6 Color-coding of the CAMLOG® Implant System according to Implant </p> <p>Please refer to the reverse.</p>
<p>7 Further Documentation</p> <p>For further information, see the current CAMLOG product catalog and the following documents:</p> <ul style="list-style-type: none">Preparation Instructions for the CAMLOG® Implant System, Art. No. 8000.0032 Instruction Manual for CAMLOG® Prosthetic Components, Art. No. J8000.0041 Instruction Manual for CAMLOG® Implant System, Art. No. J8000.0040 CAMLOG Compendium 1 Surgery, Art. No. J8000.0027 CAMLOG Compendium 2 Prosthetics, Art. No. J8000.0028 <p>The documents are available from local CAMLOG representatives or at www.camlog.com or www.camlog.de.</p>
<p>8 Contact</p> <p>Please refer to the reverse.</p>

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